

VDRIVE

Data Reduction and Interactive Visualization softwarE for Event
Mode Neutron Diffraction

Dr. Ke An

Spallation Neutron Source
Oak Ridge National Laboratory

DOCUMENT AVAILABILITY

Reports produced after January 1, 1996, are generally available free via the U.S. Department of Energy (DOE) Information Bridge:

Web site: <http://www.osti.gov/bridge>

Reports produced before January 1, 1996, may be purchased by members of the public from the following source:

National Technical Information Service

5285 Port Royal Road

Springfield, VA 22161

Telephone: 703-605-6000 (1-800-553-6847)

TDD: 703-487-4639

Fax: 703-605-6900

E-mail: info@ntis.fedworld.gov

Web site: <http://www.ntis.gov/support/ordernowabout.htm>

Reports are available to DOE employees, DOE contractors, Energy Technology Data Exchange (ETDE) representatives, and International Nuclear Information System (INIS) representatives from the following source:

Office of Scientific and Technical Information

P.O. Box 62

Oak Ridge, TN 37831

Telephone: 865-576-8401

Fax: 865-576-5728

E-mail: reports@adonis.osti.gov

Web site: <http://www.osti.gov/contact.html>

This report was prepared as an account of work sponsored by an agency of the United States Government. Neither the United States government nor any agency thereof, nor any of their employees, makes any warranty, express or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights. Reference herein to any specific commercial product, process, or service by trade name, trademark, manufacturer, or otherwise, does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States Government or any agency thereof. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States Government or any agency thereof.

VDRIVE-DATA REDUCTION AND INTERACTIVE VISUALIZATION SOFTWARE FOR EVENT MODE NEUTRON
DIFFRACTION

Dr. Ke An

SPALLATION NEUTRON SOURCE
OAK RIDGE NATIONAL LABORATORY

© 2007-2023

Prepared by
OAK RIDGE NATIONAL LABORATORY
P.O. Box 2008
Oak Ridge, Tennessee 37831-6285
managed by
UT-Battelle, LLC
for the
U.S. DEPARTMENT OF ENERGY
under contract DE-AC05-00OR22725

Preface

VDRIVE (VULCAN Data Reduction and Interactive Visualization software) is a software package that aims at time-event neutron time of flight data reduction, visualization and analysis for engineering diffractometer – the VULCAN instrument; it includes sub programs for different data reduction purposes, sequential data analysis functions, synchronization and visualization of neutron and sample environment data, and outputs of GSAS or Fullprof data for further professional data analysis. In general, it is available for event mode data reduction and analysis for TOF neutron diffractometers at spallation sources.

For the newest version of the document, it can be downloaded from the SNS data analysis computer: Applications/data analysis/VDRIVE Manual

Contact:

Dr. Ke An

kean@ornl.gov


Oak Ridge National Laboratory

Bldg 8600, MS 6475

Oak Ridge, TN 37831

Citation: when you publish your results, please cite the use of VDRIVE as:

Ke An, VDRIVE- Data Reduction and Interactive Visualization Software for Event Mode Neutron Diffraction, ORNL Report, Oak Ridge National Laboratory, ORNL-TM-2012-621, 2012.



Contents

Preface	i
1 Flowchart and cheat sheet	1
2 Some Examples	3
1.1 Generating your own instrument parameter files.	3
1.2 Single peak fit of residual strain mapping data/ in-situ step loading/heating data	4
1.3 Rietveld refinement of residual strain mapping data/ in-situ step loading/heating data	4
1.4 Single peak fit of in-situ continuous loading/heating data.....	5
1.5 Rietveld refinement of in-situ continuous loading/heating data	6
1.6 View the detector contour:.....	6
3 Computer access	6
4 VULCAN Data	6
4.1 SNS or VULCAN data folder structure	6
4.2 Instrument files location	11
4.3 Download/upload your data	11
5 Data reduction, visualization and analysis with VDRIVE	12
5.1 Load VDRIVE in IDL	12
5.2 Data Reduction and Analysis Commands	12
5.2.1 VDRIVECHOP (Chop)	12
5.2.2 VDRIVESUMGSAS (VSUM)	13
5.2.3 VDRIVEMERGE (Merge)	14
5.2.4 VDRIVEBIN (Vbin)	15
5.2.5 VDRIVEVIEW (View).....	15
5.2.6 VDRIVEFIT (VFIT).....	17
5.2.7 VPEAK.....	19
5.2.8 VDRIVEPoleP (Pole)	20
6 Data Analysis with GSAS	21
6.1 Instrument files preparation	21
6.1.1 VULCAN instrument parameter file for GSAS.....	21
6.1.2 Name conventions and the locations of instrument files	22
6.1.3 VDRIVECALI (Cali).....	22
6.1.4 VDRIVEPRM (VPRM)	23
6.2 VDRIVESPF(VSPF).....	23
6.3 VDRIVEGSAS(GSAS).....	25
Acknowledgement	27
References	27

List of Figures

Figure 1 Flowchart of VDRIVE.VDRIVE cheat sheet.....	1
Figure 3 Raw histogram plots of GSAS pattern in the two banks.....	16

Figure 3 2D contour and 3D surface plots generated by VDRIVEVIEW.....	16
Figure 5 Plot of single peak fit in one pattern by VDRIVEFIT.....	17
Figure 6 The fitting results from the VDRIVEFIT.....	18
Figure 7 Example of the VULCAN GSAS prm file.	21

1 Flowchart and cheat sheet




Figure 1 Simple flowchart of VDRIVE(X).

Note from the author:

The use of VDRIVE(X) commands in idl environment (think about python, or matlab etc.) is simple. It follows a uniform and simple grammar, i.e. command name, proposal number (IPTs), start run (runs), end run (rune), and plus some keywords if the default settings need to change. The VULCAN data storage structure is well defined on the analysis cluster, and by using the commands above VDRIVE(X) mainly requires user to provide measurements 'location' to be reduced or analyzed.

The VDRIVE commands evolve with the instrument 's upgrade including detectors, DAQ, new sample environments, and new measurement schemes, etc. New versions are introduced with major instrument upgrade. Some of the old keywords may obsolete without timely update in the manual. Let instrument team know if you have encountered any problem.

VDRIVE(X) cheat sheet

Command line mode (requires IDL license on SNS server, if not available send request to linuxsupport@sns.gov and ask for releasing taken license):

Log on SNS linux machine (<http://analysis.sns.gov>) or local machine in the hutch with your XCAMS account (The one with which you submit proposal). Right click mouse and open a **Terminal (command line)**, type **idl**, then type **@VDRIVEX** (not case sensitive for these commands anymore) for **data collected since 2021 Feb**. For data collected between 2018-2020, use **@VDRIVEbeta** add 'en' to bin and chop commands, e.g. 'vbinen'. For older legacy data, use **@VDRIVE**. The different variants are due to different detector types and upgrade.

If your data is not automatically reduced or you want to rebin the data, use the bin command below. There is often a need of runv which is a pointer to the set of calibration files including the .prm and vanadium file under the 'Instrument' folder. Refer to instrument team for the instrument file or create one using CALI command.

Note: the software is living one because of instrument upgrades, and some of the keywords are obsoleted without being updated the manual, if you encounter any problem, let instrument team know.

Bin event data to GSAS histogram files if not binned before: (short name: **vbin**)
vbin, ipts=#, runs=#, rune=# [, choprun=#]

Chop, synchronize, and bin continuously collected data in seconds: (short name: **chop**)
chop, ipts=#, runs=#, dbin=#

View one GSAS gda data pattern after binning as histogram data: (short name: **view**)
view, ipts=#, runs=# [, chopruns=#, runv =#, pc=1]

View sequential data in 2D contour and 3D surface: (short name: **view**)
view, ipts= #, runs=#, rune=# [, chopruns=#] [, minv=#.#, maxv=#.#, runv=#, norm=1 , pc=1]

Conduct GSAS single peak fit: (short name: **vspf**)
vspf, ipts=#, runs=#, rune=#, runv=# [, chopruns=#] [,runr=#] [,pc=1] [,norm=1] [, updatep=1] [, autofix=1] [, autopeak=1, npeaks=#]

*Note: peak id file should be saved as **peak.txt** as the default in the **binned_data** folder, otherwise provide the full path.*

Conduct GSAS Rietveld refinement: (short name: **gsas**)
gsas, ipts=#, runs=#, rune=# [, chopruns=#], runm=#

Generate instrument parameter and instrument spectrum files (by your local contact):

vbin, runs=#, ipts=#, tag='V'
vpeak, ipts=#, runv=#
vprm, ipts=#, runv=#, freq= 30 (60 or 20)

Generate instrument calibration files: (short name: **cali**)

cali, ipts=#, runp=#, runv=#, tag='Si', freq=20

Handy commands:

Openrecord, IPTS=12345	Open the record file
Openshare, IPTS=12345	Open the share folder
Closeplot, IPTS=12345	Close all plots
OpenMeanlog, IPTS=12345, chopruns=1111	Open the chopped log
Openlog, IPTS=12345, runs=1000	Open the collected log
Setipts, IPTS=12345	Set the working IPTS
CurrentIPTS	Get the working IPTS
IPTSprm, IPTS=12345, runv=2000	Mark the runv in IPTS sys.

VDRIVEGUI mode is not updated or supported any more. (no IDL license required, but highly recommend the command console mode):

Application/data analysis/VDRIVEGUI, or open a Terminal and type **VDRIVEGUI**.

In **VDRIVEGUI**, select **VDRIVE** functions in **VDRIVE SUB** tab to process data using the keywords above to process data.

2 Some Examples

In the examples below, the IPTS number is 1000 and the instrument parameter file is assumed to be 2000 (check instrument staff for proper instrument files for you), so make sure you have 2000-s.prm and 2000-s.gda in your "Instrument" or gda data folder.

1.1 Generating your own instrument parameter files.

Suppose a fresh calibration measurement has been done for your IPTS, you will have the 20Hz calibration powder Si run for diffraction constants, and V run for instrument spectrum. You can use a simple command to generate the instrument parameter file by the command below:

CALI, IPTS=1000, RUNP=1000, RUNV=1001

If only a change of instrument spectrum is needed due to chopper setting (e.g. 30Hz), a command below can be used to generate the specific prm file from last calibration for your IPTS.

First run **VBIN** to bin the V file:

VBIN, IPTS=1000, RUNV=1001, tag='V'

Then generate the prm file:

VPRM, IPTS=1000, RUNV=1001, freq=30

1.2 Single peak fit of residual strain mapping data/ in-situ **step** loading/heating data

Look up the runs AutoRecord.txt in the shared folder and find the correct run numbers for interest of data analysis, below suppose an IPTS 1000, with valid data from 3000 (runs, 's' means 'start') to 4000 (rune, 'e' means 'end') (missing run numbers in between is allowed), and take the 3000 as the stress free one (it can be any one in the data folder), so the data processing would be (data are now auto reduced for single runs, so you can skip the VBIN below unless you want to rebin):

VBIN, IPTS=1000, RUNS=3000, RUNE=4000

Use **VIEW** to look up the peaks you want to refine:

VIEW, IPTS=1000, RUNS=3000, RUNV=2000

Use **VIEW** to visualize the data in contour:

VIEW, IPTS=1000, RUNS=3000, RUNE=3200, RUNV=2000, NORM=1

(runv 2000 is the instrument parameter file number under the instrument folder)

Create a peak.txt file in the 'binned_data' folder (see VSPF command section for details), refine the peaks:

VSPF, IPTS=1000, RUNS=3000, RUNE=4000 [, RUNR=1] , RUNV=2000 [,Autofix=1]

Automatically fit a number of peaks by using the autopeak keyword:

VSPF, IPTS=1000, RUNS=3000, RUNE=4000 [,RUNR=1], RUNV=2000, Autopeak=1, npeaks=5

Take a look at the fitting qualities from the PDF output file to see if you need to adjust the peak positions in the peak.txt file, or use **UPDATEP=1** in **VSPF** command to update the peak position from previous fit.

1.3 Rietveld refinement of residual strain mapping data/ in-situ **step** loading/heating data

Look up the runs summary in the shared folder and find the correct run numbers for the data analysis interest, suppose you have valid run from 3000 to 4000, (allow missing run

numbers in between), and you want to take runm=9999 (yes, it can be any number even not in the existing run list) as the source one (can be any existing one in the folder, but try to not use the number within your intended refinement list because it will be overwritten) and you have two phases in the data, so the data processing would be:

VBIN, IPTS=1000, RUNS=3000, RUNE=4000

Follow the VDRIVEGSAS command details to prepare 9999_1.EXP
Perform the batch mode data refinement

GSAS, IPTS=1000, RUNS=3000, RUNE=4000, RUNM=9999, NPHASE=2

Look at the fitting qualities in the PDF, you may need to do it in segments for different setup data due to complication of phase transitions.

1.4 Single peak fit of in-situ **continuous** loading/heating data

Look up the runs summary in the shared folder and find the correct run numbers for the data analysis interest, suppose you have valid run 3000, chop the data first:

CHOP, IPTS=1000, RUNS=3000, DBIN=60, BIN=1

Suppose the data is chopped to 100 sub runs.

You can rebin the data

VBIN, IPTS=1000, CHOPRUN=3000, RUNS=1, RUNE=100

Use VIEW to look up the peaks you want to refine:

VIEW, IPTS=1000, CHOPRUN=3000, RUNS=1, RUNV=2000

Create the peak.txt file in the gda folder (see VSPF section for details), refine the peaks:

**VSPF, IPTS=1000, CHOPRUN=3000, RUNS=1, RUNE=100 [,RUNR=1], RUNV=2000
,[Autofix=1]**

Look at the fitting qualities in the PDF, you may need to adjust the peak positions in the peak.txt file, or use **UPDATEP=1** in **VSPF** command.

Automatically fit a number of peaks by using the autopeak keyword:

**VSPF, IPTS=1000, CHOPRUN=3000, RUNS=1, RUNE=100, RUNR=1, RUNV=2000,
Autopeak=1, npeaks=5**

1.5 Rietveld refinement of in-situ **continuous** loading/heating data

Look up the runs summary in the shared folder and find the correct run numbers for the data analysis interest, and you have two phases in the data, suppose you have valid run 3000, chop the data first:

CHOP, IPTS=1000, RUNS=3000, DBIN=60, BIN=1

Suppose the data is chopped to 100 sub runs.

You can rebin the data

VBIN, IPTS=1000, CHOPRUN=3000, RUNS=1, RUNE=100

Follow the VDRIVEGSAS command details to prepare 1_1.EXP

Perform the batch mode data refinement:

GSAS, IPTS=1000, CHOPRUNS=3000, RUNS=1, RUNE=100, RUNM=1, NPHASE=2

Look at the fitting qualities in the PDF, you may need to do it in segments for different setup data due to complication of phase transitions.

1.6 View the detector contour:

Pixel, IPTS=1000, RUNS=1, RUNV=2000

3 Computer access

What you need:

1 XCAMS account, the same one that you use in the IPTS system. If you do not have one follow this link to get one. (<https://user.ornl.gov>)

2 Once you have an XCAMS account, you can go to <http://analysis.sns.gov> and follow the instruction for the remote window access.

4 VULCAN Data

4.1 SNS or VULCAN data folder structure

Analyzed data are stored in a shared folder that all team members have the access. The local machine data process is disabled due to dependency on large raw neutron event and meta data.

</SNS/VULCAN/IPTS-1234/shared/>

Under this folder you may find standard folders like below

autoreduce /	auto reduced files from SNS data reduction routine. VULCAN user is recommended to use data in the <code>binned_data</code> folder below.
binned_data/	GSAS file folder by created by auto reduction or VDRIVEBIN . Some analysis results may also be in the folder along with synchronized sample environment data. Results from VDRIVESPF and VDRIVEGSAS may also store in this folder.
Pictures/Videos	Photos/Videos captured by the camera computer.
logs/	Sample environment files (loadframe, furnace and others).
Instrument/	Instrument parameter file (<code>*.prm</code>) and vanadium file (<code>*-s.gda</code>).
ChoppedData/	Chopped raw data by VDRIVECHOP . Raw data need to be reduced by running VDRIVEBIN , or selecting <code>bin=1</code> when executing VDRIVECHOP , copies of synchronized sample environment files are also in this folder.
SPF_data/	Peak fitting data created by VDRIVEFIT and VDRIVESPF .
GSAS_Results/	GSAS refinement results created by VDRIVEGSAS .
Summed_GDA /	GSAS data merged over multiple runs created by VDRIVESUMGSAS .
ViewOut/	Folder for output files when using “ <code>vout=1</code> ” keyword in VIEW .
Auto-Peak/	Created folder for storing peak files created by using “ <code>autopeak</code> ” keyword in VDRIVESPF .
AutoRecord.txt	All Record file generated by auto reduction.
AutoRecordAlign.txt	Align Data (when title has prefix of ‘Align:’) Record file generated by auto reduction.
AutoRecordData.txt	Data Record file generated by auto reduction.

Result files created by VDRIVE

[/Instrument/](#)

Under this folder you may find standard instrument parameter file pairs as below

[123456-s.gda](#) The smoothen Vanadium file, you may have more than one if different instrument optics configurations are used for your IPTS. The file can be used by specifying the run number as keyword like “`runv=123456`” as needed in different VDRIVE commands.

Vulcan-123456-s.prm The corresponding instrument parameters file including the calibration of the instrument. The file can be used by specifying the run number as keyword like “runv=123456” as needed in different VDRIVE commands.

note.txt The file contains basic information of the instrument parameter files under this folder, such as chopper or guide settings.

/SPF_data/

Under this folder you may find files created by VSPF as below

VDriveSPF-1234-1-100-bk1.pdf The snapshot of the fitting plot for each peak position fitted in VSPF. (*Runs=1 and rune =100, in chopruns=1234*)

VDriveSPF-1234-1-100-bk1.txt The output file containing fitting parameters for VSPF including a copy of the synchronized sample environment data. For example, lattice d spacings in Angstrom or TOF, integrated intensity **I** (w/o normalization to proton charge **/P**, and beam spectrum **/V**), background (**B**), **FWHM**, and its corresponding parameters, lattice strain based on either the very first runnumber or specified run number by using “runr” keyword.

The file name can be in the format below by using collected run numbers when chopping data is not required:

VDriveSPF-1000-1999-bk1.txt

/GSAS_Results/

Under this folder you may find files created by VDRIVEGSAS as below

VDriveGSAS-1234-1-100-bk1.pdf The snapshot of the fitting plot for each Rietveld refinement in **VDRIVEGSAS**. (*Runs=1 and rune =100, in chopruns=1234*)

VDriveGSAS-1234-1-100-atom.txt The output file containing atom fitting parameters by GSAS (the file with “_log” includes a copy of the synchronized sample environment data) given they are opted to be refined.

VDriveGSAS-1234-1-100-profile.txt The output file containing profile fitting parameters by GSAS (the file with “_log” includes a copy of the synchronized sample environment data) given they are opted to be refined. e.g. the weight fraction,

scale factor, lattice parameters etc. if they are opted to be refined.

The file name can be in the format below by using collected run numbers when chopping data is not required:

VDriveGSAS-1000-1999-bk1.txt

Note: The meaning of the fitted parameters can be found in the GSAS manual.

/Summed_GDA/

12345.gda Summed GSAS data merged by **VDRIVESUMGSAS**.

/ViewOut/

Under this folder you may find files created by VIEW when “vout=1” is selected.

1-123456.txt Two-column data of d spacing and Intensity as seen in the VIEW plot for bank 1, first number changes to according bank.

1-totall.txt Two-column data of runnumber and total summed Intensity as seen in the VIEW plot for bank 1, first number changes to according bank.

Note: If chopped data is processed, look for those data under the folder name of chopped run number.

/ChoppedData/1234

Under this folder you may find files created by VDRIVEChop[en].

1.dat One-column intensity data to be read by VBIN for first chopped run number. Number changes accordingly with sliced run.

choppeddatatofinfo.dat One-column tof data to be read by VBIN for all chopped run.

1234sampleenv_chopped_mean.txt Averaged synchronized sample environment data over the “dbin” time in VDRVEchop command.

1234sampleenv_chopped_start.txt Starting values of synchronized sample environment data over the “dbin” time in VDRVEchop command.

1234sampleenv_chopped_end.txt Ending values of synchronized sample environment data over the “dbin” time in VDRVEchop command.

IPTS-1000-MTSLoadFrame-1234-header.txt The header information of the sample environment data by VDRVECHOP command.

Note: If the chopped run number is 1234, for the chopped sample environment data, a copy is created in the binned_data folder/1234/.

/binned_data/1234

Under this folder you may find files created by VDRIVEChop[en].

1.gda Gsas data of all the banks. Number changes accordingly with sliced/collected run.

1234sampleenv_chopped_mean.txt Copy if averaged synchronized sample environment data over the “dbin” time in VDRVECHOP command.

1234sampleenv_chopped_start.txt Copy of starting values of synchronized sample environment data over the “dbin” time in VDRVEchop command.

1234sampleenv_chopped_end.txt Copy of ending values of synchronized sample environment data over the “dbin” time in VDRVEchop command.

IPTS-1000-MTSLoadFrame-1234-header.txt Copy the header information of the sample environment data by VDRVECHOP command.

Note: If the chopped run number is 1234, the files are found in binned_data folder/1234/.

GSAS files are in ASCII format with three columns for each bank, and both banks data are saved in one file with data structure header as a separator. The three columns are, time of flight in ms (Tof), intensity (see the file header for normalization to vanadium), square root of intensity.

The relationship between time of flight and d spacing is:

$$\text{Tof} = d * \text{DIFC} + d^2 * \text{DIFA} + \text{ZERO}$$

DIFC, DIFA and Zero are instrument parameters which can be found in the instrument parameter file (*.prm). (see instrument parameter file section for details).

Sample environment logs are in ASCII format with multiple columns. It has by default timestamp, time elapsed, and variables (loadframe signals, or temperatures).

Chopped sample environment files are synchronized with neutron timestamps. It has by default chopped run number, proton charge (for normalization by using `pcsendv=1` keyword), and synchronized columns from sample environment files.

4.2 Instrument files location

Processed instrument files and vanadium files are also in the VULCAN share folder.

/SNS/VULCAN/shared/Calibrationsfiles/Instrument/
Standards/
PRM/
Template/

4.3 Download/upload your data

Use an SFTP client to transfer data between your own computer and the data on server. FileZilla is a free one that works across different platforms. The server is analysis.sns.gov and the port for SFTP is 22.

5 Data reduction, visualization and analysis with VDRIVE(X)

5.1 Load VDRIVE(X) in IDL

VDRIVE(X) is based on the IDL.

OS	License	OPEN	Type first	Then type
Linux	Yes	Terminal	idl or idlde (case sensitive)	@VDRIVEX or vdrivex
Linux	No	Terminal	VDRIVEGUI (stop updating)	
Linux	No	Application/data analysis	VDRIVEGUI (stop updating)	

5.2 Data Reduction and Analysis Commands

5.2.1 VDRIVECHOP (Chop)

Purpose: Chop and bin continuously measured neutron data in time sequence under changing sample environment conditions.

Common use:

CHOP, IPTS=1000, RUNS=2000, dbin=60, bin=1

where, **dbin** is the chop step size in seconds; **bin=1**, for binning data to GSAS file after slicing the data in time. GSAS data are stored at **/SNS/VULCAN/IPTS-1000/shared/binned_data/2000/** along with the chopped sample environment files **2000sampleenv_chopped_start(mean or end).txt**.

Alternate for **loadframe=1: furnace=1**, or **generic=1**, when using VULCAN standard sample environment DAQ for the furnaces or others. For a customized sample environment file name, use **SampleEnv='your sample file name.txt'** (the customized sample environment file is stored in **/SNS/VULCAN/IPTS-1000/shared/logs**).

If no sample environment is chosen or **justchop=1** keyword is selected, no sample environment data synchronization will be executed.

Other uses:

To chop with customized time segments:

CHOP, IPTS=1000, RUNS=2000, pickdata="/SNS/VULCAN/IPTS-1000/shared/picktime.txt", bin=1, pulsetime=1

where **pickdata** is the file name containing several selected time segments of the neutron data with the format of below (separated by tab and unit is seconds and use pulse time):

1.00 10.00
20.01 30.00
40.01 50.00

...

Additional keywords:

Focus_EW=0	Bin data over each detector module.
dt=###.##	Time (s) between each chopped run, by default it equals dbin . If less than dbin , each run will have overlapped neutron events.
t0=###.##	Time (s) offset from the start of neutron event files.
tend=###.##	Time (s) offset from the end of neutron event files.
Ndataset=###	Number of chopped data set. Effective when the automatically calculated data sets are more than Ndataset .
Ndbin=###	Number of time bins per cycle for stroboscopy.
Ncycle=###	Number of cycles to perform stroboscopy.
Strobe=1	Set stroboscopy on.
GSAS=1	Plot GSAS file after binning.
Onelog=1	Read sample environment parameters from one log file.
Loadframe=1	Read sample environment parameters from MTS Loadframe log file.
Generic=1	Read sample environment parameters from Generic log file e.g. analog inputs, voltage, current, exclusive with loadframe=1.
Accu=1	Accumulate data over chopped runs.
PulseTime=1	For PickData , when PulseTime is used in the picktime.txt file.
AppRun=####	Append chopped data to previously chopped data in choprun=#### . Will not work with Connect=1 .
Connect=1	Chop more than one runs and connect data after RUNS and data are saved to RUNS's chopped data folder. Will not work with AppRun .

Other **VDRIVEBIN** keywords when coupled with **bin=1**.

5.2.2 VDRIVESUMGSAS (VSUM)

Purpose: Combine GSAS data. GSAS data are combined from the runs.

Common use:

VSUM, IPTS=1000, RUNFILE="/SNS/VULCAN/IPTS-1000/shared/runfile.txt"

VSUM, IPTS=1000, RUNLIST=[1234,1235]

VSUM, IPTS=1000, RUNS=1234,RUNE=1235

The summed GSAS data are saved to [/SNS/VULCAN/IPTS-1000/shared/Summed_GDA/](#)

Example of the tab delimited runfile.txt:

```
-----  
1001 1002 1003 1004  
1005 1006 1007  
1008 1009 1010  
...  
-----
```

Additional keywords:

Runv=1111 Normalize summed data over Vanadium.
Trans=1 Sum data to the first row runs.
Threeframes=1 Sum data measured by different chopper settings.

5.2.3 VDRIVEMERGE (Merge)

Purpose: Combine collected event data (*not work with new nED data anymore*). Data are combined from the runs of rest columns to the runs of the first column in the runfile.txt.

Common use:

**MERGE, IPTS=1000, RUNFILE="/SNS/VULCAN/IPTS-1000/shared/runfile.txt",
CHOPRUN=2**

The combined data are saved to [/SNS/VULCAN/IPTS-1000/shared/chopped_data/2/](#)

To bin the data combined by VDRIVEMERGE:

VDRIVEBIN, IPTS=1000, CHOPRUN=2

GSAS files are stored in [/SNS/VULCAN/IPTS-1000/shared/binned_data/2/](#)

Example of the tab delimited runfile.txt:

```
-----  
1001 1002 1003 1004  
1005 1006 1007  
1008 1009 1010  
...  
-----
```

Additional keywords:

NONE

5.2.4 VDRIVEBIN (Vbin)

Purpose: Bin collected data or chopped data by **VDRIVECHOP** to GSAS files.

Common use:

For typical mapping experiments or single run:

VBIN, IPTS=1000, RUNS=2000, RUNE=2099

GSAS files are stored in [/SNS/VULCAN/IPTS-1000/shared/binned_data/](#)

*For chopped files created by **VDRIVECHOP**:*

VBIN, IPTS=1000, CHOPRUN=2000, RUNS=1, RUNE=100

GSAS files are stored in [/SNS/VULCAN/IPTS-1000/shared/binned_data/2000/](#)

Additional keywords:

BINW=0.005

The logarithm bin step size of TOF and its default is 0.001.
For 7-module data binning.

FOCUS_EW=0

RUNV=5000

Normalize data over smoothed vanadium file in
[/Instrument](#) folder, and used for **VDRIVEFit** when
normalized intensity is an output option, for **VDRIVESPF**,
NORMALIZE=0 should be taken.

Runb=5001

Correct data by background file, use with **NORMALIZE=1**

IParm='prmfile.prm'

Replace the **Vulcan.prm** in comment line with the
customized **prmfile.prm** string in GSAS files.

FullProf=1

Output FullProf files.

NoGSAS=1

Omit GSAS files.

PlotFlag=1

Plot histogram after binning.

Merge_EW=1

Bin east and west banks data into bank 1 histogram.

NoMask=1

Bin bad pixels too.

Tag='Si'

Bin a Si calibration powder. If 'CeO2' is chosen, bin a CeO2
calibration powder. If 'V' is chosen, bin a vanadium data.

Mytofbmax=1111

max tof bin in us

Mytofbmin=1111

min tof bin in us

5.2.5 VDRIVEVIEW (View)

Purpose: Visualize data in 2D, 3D plots.

Common use:

To view one GSAS pattern:

VIEW, IPTS=1000, RUNS=1




Figure 2 Raw histogram plots of GSAS pattern in the two banks.

To view sequential GSAS patterns in 2D and 3D:

VIEW, IPTS=1000, RUNS=1, RUNE=50

To view sequential GSAS patterns of chopped data in 2D and 3D:

VIEW, IPTS=1000, CHOPRUN=2000, RUNS=1, RUNE=50




Figure 3 2D contour and 3D surface plots generated by VDRIVEVIEW.

Additional keywords:

RUNV=5000

Normalize GSAS pattern with instrument spectrum by the vanadium pattern.

TOTALCOUNTS=1

Normalize with total counts of each bank.

PCSENV=1

For chopped data, normalize chopped data with proton charge (beam charge).

MinV=0.7, MaxV=2.7

Define the display range in d space.

Position=1, SENV=5

Show the 6th variable (starting from 0) as the vertical axis in the contour and surface plots. Can be any Nth column in the chopped sample environment files.

Norm=1	Normalize to proton charge of single run measurement (differs from PCSENV for chopped runs).
Tof=1	Show x in time of flight.
Q=1	Show x in Q.
PlotGSAS=1	Show individual GSAS pattern.
ClearPlot=1	Reset plotted windows by previous execution of VDRIVEVIEW .
SameScale=1	Show two banks data in same scale.
Merge=1	Show two banks data in one window.
Log=1	Plot intensity in logarithm.
Summed=1	View summed GSAS data in Summed_GDA folder.

5.2.6 VDRIVEFIT (VFIT)

Purpose: Gaussian function single peak fit and results visualization.

Common use:

Note: use *VDRIVESPF* whenever possible.

For one GSAS pattern fit:

VFIT, IPTS=1000, RUNS=1, listd=[2.60,2.45,1.89,1.6,1.47], width=[0.035, 0.03, 0.03, 0.03, 0.03], plot=1

where, **listd** is the list of initial guess of the peak position, **width** is the data range of the peak in d-space for each peak, can be one value for all peaks, or an array of values for each peak;




Figure 4 Plot of single peak fit in one pattern by VDRIVEFIT.

For sequential GSAS patterns fit:

VFIT, IPTS=1000, RUNS=1, RUNE=50, RUNR=1, listd=[2.04, 1.76, 1.243, 1.06], width=[0.03,0.03,0.03,0.02], UpdateP=1

or for chopped GSAS patterns fit

VFIT, IPTS=1000, ChopRun=2000, RUNS=1, RUNE=50, RUNR=1, listd=[2.04, 1.76, 1.243, 1.06], width=[0.03,0.03,0.03,0.02], UpdateP=1

where **RUNR** is the strain reference run, otherwise, d-space will be the only output; **UpdateP=1**, peak center of previous run will be used as the guess of current run.




Figure 5 The fitting results from the VDRIVEFIT.

This program can search peaks based on statics:

VFIT, IPTS=1000, RUNS=1, RUNE=50, RUNR=1, npeaks=8, UpdateP=1

Additional keywords:

Bterm=1	Background parameter. Default is 1 for flat background, 2 is for linear background.
Bank=[1,2]	Choose specific bank data to fit.
Norm=1	Normalize the integrated intensity to total counts.
Voigt=1	Peak shape function option. Default peak shape is Gaussian, if Voigt=1 , VDRIVEFIT will use the asymmetric Voigt function (Note: Voigt function has 2 more parameters than Gaussian function, so if the peak intensity is low, the fit may not converge).
NoErrorBar=1	Omit the error bars in plots.
Position=1, SENV=5	Show the 6th variable (starting from 0) as the vertical axis in the contour and surface plots. Can be any Nth column in the chopped sample environment files.
Log=1	Take logarithm of the intensity and then perform the fit.
Pcsenv=1	Normalize the intensity by the proton charge in the chopped sample environment file.
Sho=1	Show the fitting results in the command console.
Showbad=1	Show the bad fitting results in plots.
Npeaks=1	Set number of highest peaks (raw data) to be detected and conduct fitting by the program itself.
AutoPeak=1	Output the fit result (only with one run fitting) to a peak txt file for VSPF input.

Peak positions, intensities, peak widths and strains of each peak are stored in ASCII files **RUNS_1.txt** and **RUNS_2.txt** under **/SNS/VULCAN/IPTS-####/shared/SPF_data**. Results should be checked. Usually if there are some exotic points, it means either the statistics of these data are poor, or some initial parameters such as peak position or width should be adjusted. If peaks are too close, initial values of original position and width are important. This command is not recommended for overlapping peaks fitting.

5.2.7 VPEAK

Purpose: Process vanadium diffraction peak and data noise.

Common use:

VPEAK, IPTS=1000, RUNV=5000

Additional keyword:

Nsmooth =51	The number of points to be used in the boxcar smoothing algorithm, the bigger the smoother.
OneBank=1	When all banks' data are binned as one bank data.
Shift=1	When the chopper center is shift to large lambda aggressively.

The smoothed data is named as **####-s.gda** and located at **/SNS/VULCAN/IPTS-1000/shared/Instrument** as well as a copy in the VULCAN shared folder **/SNS/VULCAN/shared/Calibrationsfiles/Instrument/Standards/**.

5.2.8 VDRIVEPoleP (Pole)

Purpose: Process single peak data texture measurement and output pole figure data.

Common use:

VDRIVEPoleP, IPTS=1000, RUNS=1, RUNE=120, bina=1, focus_ew=0

Additional keyword:

Involume=1	When the sample is a thin and long bar measured. If the sample fully embraces the 5x5x5mm ³ gauge volume, no need to use this keyword.
Std4590=1	When the measurement uses standard rotation of Omega or VROT from 45 to 90 degrees.
rotation=value of angle	Specify a value of angle in degree to rotate the plot, default is -90.
bina=1	Bin the data with overlaps in pole/Q space, and plot will be generated if the binning is successful.

A detailed instruction of generating texture pole figure can be found in Appendix.

6 Data Analysis with GSAS

6.1 Instrument files preparation

6.1.1 VULCAN instrument parameter file for GSAS

```
1 30hz HI customized 3x15x2mm slit
2 ID Si Standard Run 6886 Vanadium 6887, 30hz HI center=2.0, band=2.88
3 INS BANK 2
4 INS FPATH1 43.754
5 INS HTYPE PNTR
6 INS NSPEC 6
7 INS 1 ICONS 16364.280 -0.2200 0.0000 0.000 0 0.000
8 INS 1 IRAD 0
9 INS 1BNKPAR 1.5000 -90.00 0.00 .00000 .3000 1 1
10 INS 1I HEAD No title
11 INS 1I ITYP 10 6.5000 42.0000 42432
12 INS 1INAME Vulcan
13 INS 1PRCF1 3 21 0.00050
14 INS 1PRCF11 1.0 0.0601700 0.00203650 0.000000E+000
15 INS 1PRCF12 612.000000 0.000000E+000 0.000000E+000 4.040000
16 INS 1PRCF13 0.000000E+000 0.000000E+000 0.000000E+000 0.000000E+000
17 INS 1PRCF14 0.000000E+000 0.000000E+000 0.000000E+000 0.000000E+000
18 INS 1PRCF15 0.000000E+000 0.000000E+000 0.000000E+000 0.000000E+000
19 INS 1PRCF16 0.000000E+000
20 INS 2 ICONS 16368.300 -0.870 0.0000 0.000 0 0.000
21 INS 2 IRAD 0
22 INS 2BNKPAR 1.5000 90.00 0.00 .00000 .3000 1 1
23 INS 2I HEAD No title
24 INS 2I ITYP 10 6.5000 42.0000 40348
25 INS 2INAME Vulcan
26 INS 2PRCF1 3 21 0.00050
27 INS 2PRCF11 1.0 0.0615814 0.00873451 0.000000E+000
28 INS 2PRCF12 612.100000 0.000000E+000 0.000000E+000 5.100000
29 INS 2PRCF13 0.000000E+000 0.000000E+000 0.000000E+000 0.000000E+000
30 INS 2PRCF14 0.000000E+000 0.000000E+000 0.000000E+000 0.000000E+000
31 INS 2PRCF15 0.000000E+000 0.000000E+000 0.000000E+000 0.000000E+000
32 INS 2PRCF16 0.000000E+000
33 INS 1 MFIL 6887-s.gda
34 INS 2 MFIL 6887-s.gda
35
```

Figure 6 Example of the VULCAN GSAS prm file.

Contents of instrument (*prm) file

Basically the prm file contains 2 parts:

1) Instrument spectrum for normalization

Lines with keyword “ITYP” and “MFIL” are related to instrument spectrum file by vanadium. In the ITYP line, “10” means the normalization method is using vanadium pattern directly, “6.5” and “42” are the valid data range in TOF ms. The MFIL lines specify the vanadium file that is used.

2) Instrument calibration and peak profile information

Lines with keyword “ICONS” and “PRCF” are related to the calibration information of the instrument. In the ICONS line, DIFC, DIFA and Zero are given. In the PRCF line, the

peak profile parameters are given. These values should be obtained by fitting the pattern of a standard sample e.g. Si.

More information of the instrument parameter files can be found in the GSAS manual.

If needed, instrument files can be created/updated by refining the standard data Si every cycle along with the smoothed vanadium data. Or one can look for one instrument file that has the same configuration (guide, chopper, lambda center, lambda width, and sample environment) from [/SNS/VULCAN/shared/Calibrationfiles/Instrument/PRM](#). The files are updated usually once in a cycle.

6.1.2 Name conventions and the locations of instrument files

VULCAN GSAS instrument files are named according to the vanadium run number, such as [Vulcan-5000-s.prm](#), where 5000 is the vanadium run number, and “-s” means the vanadium pattern has been smoothed by **VPEAK**.

For each user project, the instrument files, including the prm file and the smoothed vanadium data, should be stored in the IPTS folder [/SNS/VULCAN/IPTS-####/shared/Instrument](#).

The template prm files are in the folder [/SNS/VULCAN/shared/Calibrationfiles/PRM/Template](#).

Note: User may have more than one instrument prm files depending on the configurations used in the experiment.

6.1.3 VDRIVECALI (Cali)

Purpose: Generate instrument calibration files from Si and vanadium measurements.

Common use:

CALI, IPTS=4744, RUNP=12474, RUNV=12475, TAG='Si', Freq=20

This example is based on high intensity (HI) 20Hz chopper setting Si powder measurement **RUNP=12474**, and the corresponding vanadium measurement **RUNV=12735**. Instrument parameter file [Vulcan-12735-s.prm](#) will be created in [/SNS/VULCAN/shared/Calibrationfiles/Instrument/PRM/](#) and also, the template instrument file [Vulcan-template-HI.prm](#) is created or updated under [/SNS/VULCAN/shared/Calibrationfiles/Instrument/Template](#)

Additional keywords:

OneBank=1

For one bank calibration when two banks data are binned as one. The files will be saved in ['/1bk'](#) under those instrument folders.

6.1.4 VDRIVEPRM (VPRM)

Purpose: Generate user specific instrument files from calibration (template) files.

Common use:

VPRM, IPTS=1000, RUNV=5000, FREQ=20 (or 30, 60)[, ONEBANK=1]

Vulcan-5000-s.prm and 5000-s.gda will be created in /SNS/VULCAN/IPTS-1000/shared/Instrument.

6.2 VDRIVESPF(VSPF)

Purpose: Use GSAS for single peak fit including overlapping peaks.

Common use:

For typical mapping experiments:

VSPF, IPTS=1000, RUNS=1, RUNE=100, RUNR=1, RUNV=5000

For sequential GSAS files of chopped run:

VSPF, IPTS=1000, CHOPRUN=2000, RUNS=1, RUNE=100, RUNR=1, RUNV=5000, PCSENV=1

Required files: Peak ID file (**peak.txt**), instrument parameter file (***.prm**), and instrument spectrum file by vanadium (***-s.gda**). ***.prm** and ***-s.gda** should be under the instrument folder. The peak ID file is recommended to be named as **peak.txt** under the (chopped) GSAS data folder.

An example of the **peak.txt** file is below:

```
-----  
$bk  name  number      pos      range  
1    110   1        2.02692  0.0300000  
  
$1   310   1        0.906467 0.0250000  
  
$2   110   1        2.02692  0.0300000  
  
$2   200   1        1.43325  0.0300000  
  
2    211   1        1.17024  0.0200000  
  
$2   220   1        1.01346  0.0300000  
-----
```

The names of the columns are:

Bank ID, name of the peak, number of peaks, estimated peak position (in d), estimated peak range (in d). \$ sign is for comment line.

For overlapped peaks, which are too close to perform single peak fit, alternate the first peak column value for the purpose of fitting and outputting the corresponding results. An example is given below:

```

-----
$ Fit the first peak
1   peak1 2   2.02345           2.03456           0.04
$ Fit the second peak
1   peak2 2   2.03456           2.02345           0.04
-----

```

Additional keywords:

peakFile='peak.txt'	For customized peak ID file name other than the default one.
Runfile='file.txt'	Fit runs in a text file, one run per line.
UdataP=1	Take previous run as the peak parameters guess for current run.
prmFile=' prmfile.prm'	Use different prm file.
myRefine='Refine.txt'	Use a customized refinement steps in a macro file.
NoErrorBar=1	Omit the error bars in plots.
UserDataDir='/folder'	Direct the data in a specific folder.
Plotdata=0	Omit the plots.
Showbad=1	Show bad refinement results in the plots.
Pcsenv=1	Spectrum divided by proton charge, not with monitor .
Monitor=1	Spectrum divided by monitor 2, not with pcsenv .
Normalization=0	Not normalize to Vanadium file (with runv) for data already normalized by VDRIVEBIN.
AutoPeak=1	With the keyword below to automatically search for number of peaks from the order of high to low intensity and allow user to fit the data with the found peaks; not with "peakfile" .
NPeaks=#	Specify number of peaks to search in Autopeak option.
AutoFix=1	When a fitting fails due to a miss input of the initial peak position, this command enable the program to automatically adjust the input value over small step (dtweack=0.001 ,) over a total number of iterations (ntweack=200) to fit the data.

Result file list:

VDriveSPF-2000-1-50-bk1.txt	Peak positions, widths, intensities, and strains (in refinement folder as well as a copy in SPF_data folder).
VDriveSPF-2000-1-50-bk1.pdf	Fitting plots for quality check (in refinement folder as well as a copy in SPF_data folder).
VDriveSPF-2000-1-50-bk1.log	Refinement histories.

6.3 VDRIVEGSAS(GSAS)

Purpose: Use GSAS for Rietveld Refinement based on a source run.

Common use:

For typical runs:

GSAS, IPTS=1000, RUNS=1, RUNE=100, RUNM=9999, BANK=1

For GSAS files from chopped data:

GSAS, IPTS=1000, CHOPRUN=2000, RUNS=1, RUNE=100, RUNM=9999, BANK=1

Required files: Instrument parameter file (*.prm), instrument spectrum file by vanadium (*.s.gda) and a source run (RUNM) which has been refined well by GSAS are required and stored the GSAS data folder before the execution of this command. The source files will define the refinement scheme of Rietveld refinement in GSAS. If the source RUNM is 9999, VDRIVEGSAS needs GSAS EXP files 9999_1.EXP and 9999_2.EXP for bank 1 and 2, respectively.

Additional keywords:

Nphase=2	Number of phases in the GSAS data.
UserDataDir=' /folder'	Direct the data in a specific folder.
Runfile='runlistfile.txt'	Fit runs in a text file, one run per line.
Title='Title'	Title for all GSAS refinement.
Samerunm=1	Use the same source runm for all GSAS refinement
Extractlist=1	Extract data from lst file only, use after manually refining individual runs.

Result file list:

VDriveGSAS-2000-1-50-bk1.txt	Lattice parameters, strains.
VDriveGSAS-2000-1-50-bk1.pdf	Refinement plots for quality check
VDriveGSAS-2000-1-50-bk1-atom.txt	Atom occupancies if turned on.
VDriveGSAS-2000-1-50-bk1-profile.txt	Peak profile parameters.
VDriveGSAS-2000-1-50-bk1.log	Refinement histories.

Create a source run for VDRVIEGSAS

Open Terminal

Type *expgui*

Direct to your "binned_folder" then type #####_1.EXP,(##### can be any number even not exist in the current runlist, which is recommended) then click "Create" button

Type a string for the exp file title.

Click "Add Phase"

Select the way to add new phase (either by previous EXP file, or CIF file)

Point to the file containing crystal phase information

Click "Continue"

Modify the atom name properly and click "Add Atoms"

In "Powder" tab, click "Add New Histogram"

In Data file, click "Select"

In Instrument file, click "Select" (better copy the instrument files to the binned_data folder)

Follow the GSAS manual to perform a good refinement of the data.

Perform same way to generate source run for bank 2 by simply change *####_1.EXP* to *####_2.EXP*.

Acknowledgement

The computer and software support from Harley Skorpenske, Ling Yang, Jean C. Bilheux, John Quigley, Kyle Anderson, and Jason Hodges are acknowledged. The Spallation Neutron Source is a user facility at Oak Ridge National Laboratory sponsored by the Office of Basic Energy Science, the Office of Science, the Department of Energy, USA. Oak Ridge National Laboratory is managed by UT-Battelle, LLC, for the U.S. Dept. of Energy under contract DE-AC05-00OR22725.

References

1. A.C. Larson and R.B. Von Dreele, "General Structure Analysis System (GSAS)", Los Alamos National Laboratory Report LAUR 86-748 (2000).
2. IDL, <http://www.exelisvis.com/>

Appendix

Texture pole figure instruction on VULCAN

9.22.20 (for data collected before 2021 Feb)

For texture data collected after 2021 Feb, see next instruction.

1. load VDRIVE

idl, @vdrivebeta

2. Create peak.txt file for fitting of 6 modules/detectors and save in the binned_data folder

For cubic system first 4 peaks are good start. Note here: **no underscore** should be used in peak id. For FCC you can just choose a few first peaks, and fit same peaks in six banks.

peak.txt				
\$bk	Pkid	Npk	pos	range
1	111B1	1	2.34	0.03
1	200B1	1	2.02	0.03
#1	220B1	1	1.43	0.03
#1	311B1	1	1.22	0.03
#1	222B1	1	1.16	0.03
#1	400B1	1	1.01	0.03
#1	331B1	1	0.93	0.025
#1	420B1	1	0.90	0.025
2	111B2	1	2.34	0.03
2	200B2	1	2.03	0.03
#2	220B2	1	1.43	0.03
#2	311B2	1	1.22	0.03
#2	222B2	1	1.16	0.03
#2	400B2	1	1.01	0.03
#2	331B2	1	0.93	0.025
#2	420B2	1	0.90	0.025
3	111B3	1	2.34	0.03
3	200B3	1	2.02	0.03
#1	220B1	1	1.43	0.03
#1	311B1	1	1.22	0.03
#1	222B1	1	1.16	0.03
#1	400B1	1	1.01	0.03
#1	331B1	1	0.93	0.025
#1	420B1	1	0.90	0.025

We used standard rotation process to measure texture data, i.e. rotate vertical axis 45 degrees by 5 degree/step, and rotate horizontal axis by 360 degrees by 30 degree/step. The total measurements are 120 runs.

3 For a specific sample, find the first run with B1=0, B2=-90

187487	IPTS-25434	Sample 5	B1=0, B2=-90	b'
187488	IPTS-25434	Sample 5	B1=0, B2=-90	b'
187489	IPTS-25434	Sample 5	B1=0, B2=-90	b'

4 Find the last run with B1=45, B2=-45

187604	IPTS-25434	Sample 5	B1=45, B2=-45	i
187605	IPTS-25434	Sample 5	B1=45, B2=-45	i
187606	IPTS-25434	Sample 5	B1=45, B2=-45	i

5 Rebin the 120 data runs into individual banks by using the keyword focus_ew=0 in VBINEN

```
IDL> vdrivebinen, ipts=25434, runs=187487, rune=187606, focus_ew=0
```

6 Do the single peak fitting with VSPF as usual except that the runv needs to contain 7-bank calibration, ask instrument staff if it is not available

```
IDL> vspf, ipts=25434, runs=187487, rune=187606, runv=187154
```

Double check the fitting quality for peaks showing, and adjust the input peak values in the peak.txt file, 'update=1' is not recommended as in certain angle, for a strong textured sample, peaks would be missing.

7 Run VDRIVEPOLEP to process the data for pole figures

```
vdrivepolep, ipts=25434, runs=187487, rune=187606, involume=1, focus_ew=0  
[std4590=1]
```

The use of "involume=1" is when you have a thin and long bar measured. If the sample is large enough contains fully the 5x5x5mm³ gauge volume, or the sample is a small cube and fully buried within the gauge volume, no need to use "involume=1".

Use "std4590=1" when the measurement is conducted by standard rotation of Omega or vertical rotation axis from 45 degree to 90 degree.

This command will simply output 120 records after coordinates transformation and intensity normalization. At this point the output data as measured can be plotted in plotting software.

Below is recommended to run with "bina=1" for binning the data with overlaps in pole/Q space, and plot will be generated if the binning is successful.

vdrivepolep, ipt=25434, runs=187487, rune=187606, involume=1, focus_ew=0, bina=1 [std4590=1]

Finally, the file outputs are


pole data is saved in: /SNS/VULCAN/IPTS-25434/shared/SPF_data/VDriveSPF-187487-187606_pole.txt
 Sample Psi and Phi data is saved in: /SNS/VULCAN/IPTS-25434/shared/SPF_data/VDriveSPF-187487-187606_Psi_Phi.txt
 a copy of the pole jpg is saved in: /SNS/VULCAN/IPTS-25434/shared/SPF_data/VDriveSPF-187487-187606_pole.jpg

Name	Size	Modified
VDriveSPF-187487-187487-bk1.pdf	18.4 kB	20:34
VDriveSPF-187487-187487-bk1.txt	3.4 kB	20:34
VDriveSPF-187487-187606-bk1.pdf	905.2 kB	20:38
VDriveSPF-187487-187606-bk1.txt	225.9 kB	20:38
VDriveSPF-187487-187606_pole.jpg	538.2 kB	21:19
VDriveSPF-187487-187606_pole.txt	12.2 kB	21:19
VDriveSPF-187487-187606_Psi_Phi.txt	12.2 kB	21:19


A quick view of the pole figure of the fitted peak is shown, and the files are saved in the **SPF_data** folder *. _pole.txt is the data of *_pole.jpg file. Use *_Psi_Phi.txt to plot pole figure in other plotting software, such as OriginPro.

When using bina=1, a pole plot is generated as shown by the example below.

VULCAN Instrument layout



Schematic illustration of the bead-on-plate experimental set-up on VULCAN (top view, not to scale). The -90° and $+90^\circ$ detector banks record diffraction peaks of the (h k l) lattice planes whose normals are parallel to Q_1 and Q_2 , respectively. Strain components along these two directions are measured simultaneously. The bead-on-plate specimen is positioned on top of the sample stage and aligned at 45° from the incident beam.



Instrument layout and an example pole figure from the fitted data.

In the pole figure's polar coordinates, the radial coordinate is linked to Psi (or omega of the instrument), and the angular coordinate is linked to Phi (or HROT of the instrument motor). With relative to the sample mount starting as above, Q_1/LD is at the center of

the pole middle (pointing out), Q2/ND at 0 degree, and Q3/TD (perpendicular to the paper in the instrument plot) is at -90 degree.

Currently the plotted pole figure is rotated about -90 degree from the measurement. You can adjust it back by using "rotation=0" in the VDRIVEPOLEP command. This can also be used to slightly adjust the offset of the rotation by specifying an angle in degree. When using "bina=1", adjusting rotation with this keyword may produce error in area binning. You can tweak it to avoid it.

New Texture pole figure instruction on VULCAN (beta)

For data collected after 2021 Feb

1. load VDRIVE

idl, @vdrivex

2 For a specific sample, find the first run with B1=0, B2=-90

187487	IPTS-25434	Sample 5	B1=0, B2=-90	b'
187488	IPTS-25434	Sample 5	B1=0, B2=-90	b'
187489	IPTS-25434	Sample 5	B1=0, B2=-90	b'

3 Find the last run with B1=45, B2=-45

187604	IPTS-25434	Sample 5	B1=45, B2=-45	i
187605	IPTS-25434	Sample 5	B1=45, B2=-45	i
187606	IPTS-25434	Sample 5	B1=45, B2=-45	i

4.Holistic binning, fitting, pole data, processing and plotting.

Texture, `ipts=26807`, `runs=firstrun`, `rune=lastrun`, `runv=texturevanadiumrun`, `listd=[peakposition1, peakposition2,...]`, `width=[width1, width2,...]`, `/autopeak`

Note: number of elements in listd and width should be equal. They are initial estimate and autopeak will refine them spit out a peak data for vdrivespf. Check on the file to decide to continue the process or not.

You can then play with VDRIVEPOLEP command to change the different parameters as the prior instruction.