

PDF reduction

(no inelastic, absorption or multiple
scattering correction for sample,
correction is done for V rod)

1. In mantidplotnightly use 'LoadGSS' to load autoreduced gsas file
2. 'ConvertUnits' to momentum transfer
3. 'Scale' by Multiplying the data by 0.01
4. Plot Spectrum
5. Fit Q space from 30 – 50 using a flat background
6. Divide the spectrum by the flat background coefficient A0 so that the q data $\rightarrow 1$ using scale command. This is your arbitrarily normalized $S(Q)$
7. Now use 'PDFFourierTransform' to chose Qmin, Qmax etc. to produce $G(r)$
8. SavePDFGui to generate Gr for PDFGui which generates relevant header for neutron.

because there is no connection to the host reports.mantidproject.org
Hint: Check your connection following this link: <http://reports.mantidproject.org/api/feature>><http://reports.mantidproject.org/api/feature> Host not found: reports.mantidproject.org

Algorithms

Execute LoadGSS

Algorithms

- Arithmetic
- CorrectionFunctions
- Crystal
- DataHandling
- Diagnostics
- Diffraction
- Events
- Examples
- ILL
- Inelastic
- MDAgorithms
- Muon

Workspaces

Load ▼ Delete Group Sort ▼ Save ▼

Filter Workspaces

Workspaces

LoadGSS input dialog

Loads a GSS file such as that saved by SaveGSS. This is not a lossless process, as SaveGSS truncates some data. There is no instrument associated with the resulting workspace. 'Please Note': Due to limitations of the GSS file format, the process of going from Mantid to a GSS file and back is not perfect.

Filename

roups/pg3-cycles/2017-A-B-data/PDF/PG3_37883.gsa

Browse

OutputWorkspace

BaTiO3

UseBankIDasSpectrumNumber

?

Keep Open

Run

Close

ConvertUnits input dialog

Performs a unit change on the X values of a workspace

InputWorkspace: BaTiO3

OutputWorkspace: BaTiO3_Q

Target: MomentumTransfer

EMode: Elastic

EFixed:

☐ AlignBins

☒ ConvertFromPointData

? Keep Open ☐ Run Close

Y is multiplied by bin width

L1 = 60

Bank 0: Total flight path = 63.18 2Th

LoadGSS successful, Duration 0.08 sec

Algorithms

Execute ConvertUnits

Algorithms

- Arithmetic
- CorrectionFunctions
- Crystal
- DataHandling
- Diagnostics
- Diffraction
- Events
- Examples
- ILL
- Inelastic
- MDAlgorithms
- Muon

Workspaces

Load ▾ Delete Group So

Filter Workspaces

Workspaces

- BaTiO3

Scale input dialog

Scales an input workspace by the given factor, which can be either multiplicative or additive.

InputWorkspace	BaTiO3_Q	⌚
OutputWorkspace	BaTiO3_Q_scaled	🗑️
Factor	0.01	🗑️
Operation	Multiply	⬆️⬆️

? Keep Open ☐ Run Close

Name: ConvertFromPointData
ConvertUnits successful, Durati

Algorithms

Execute Scale

Algorithms

- ▶ Arithmetic
- ▶ CorrectionFunctions
- ▶ Crystal
- ▶ DataHandling
- ▶ Diagnostics
- ▶ Diffraction
- ▶ Events
- ▶ Examples
- ▶ ILL
- ▶ Inelastic
- ▶ MDAlgorithms
- ▶ Muon


Workspaces

Load ▼ Delete Group

Filter Workspaces

Workspaces

- ▶ BaTiO3
- ▶ BaTiO3_Q



DejaVu Sans 10 B /t U x² x₂ αβ Γ f


Fit Function (Chi-sq = 2.30391, failed)

Fit

Display

Setup

Property	Value
▼ Functions	
Type	CompositeFunction
NumDeriv	☐ False
▼ f0-FlatBackground	
Type	FlatBackground
A0	2.585833
▼ Settings	
Workspace	BaTiO3 Q scaled
Workspace Index	0
StartX	30.000000
EndX	55.000000
Output	BaTiO3 Q scaled
Minimizer	Levenberg-Marquardt
Ignore invalid ...	☐ False
Cost function	Least squares
Max iterations	500
Peak Radius	0
Plot Difference	☒ True
Plot Composite ...	☐ False
Convolve Comp...	☐ False
Show Paramete...	☐ False
Evaluate Functi...	CentrePoint



- Add peak...
- Add background...
- Add other function...
- Plot guess
- Reset range
- Clear all
- Get Parameters
- Fit
- Undo fit

Results Log

Parameters:

Name: InputWorkspace, Value: BaTi

Name: OutputWorkspace, Value: Ba

Direction: Output

Name: Factor, Value: 0.01, Default?

Name: Operation, Value: Multiply, I

Scale successful, Duration 0.00 secon

Algorithms

Execute

Scale

Algorithms

- ▶ Arithmetic
- ▶ CorrectionFunctions
- ▶ Crystal
- ▶ DataHandling
- ▶ Diagnostics
- ▶ Diffraction
- ▶ Events
- ▶ Examples
- ▶ ILL
- ▶ Inelastic
- ▶ MDAlgorithms
- ▶ Muon

Workspaces

Load

Delete

Group

Sc

Filter Workspaces

Workspaces


- ▶ BaTiO3
- ▶ BaTiO3 Q
- ▶ BaTiO3 Q scaled

Fit Function (Chi-sq = 2.30391, failed)

Fit Display Setup

Property	Value
▼ Functions	
Type	CompositeFunction
NumDeriv	☐ False
▼ Settings	
Workspace	BaTiO3 Q scaled
Workspace Index	0
StartX	30.000000
EndX	55.000000
Output	BaTiO3 Q scaled
Minimizer	Levenberg-Marquardt
Ignore invalid ...	☐ False
Cost function	Least squares
Max iterations	500
Peak Radius	0
Plot Difference	☒ True
Plot Composite ...	☐ False
Convolve Comp...	☐ False
Show Paramete...	☐ False
Evaluate Functi...	CentrePoint

BaTiO3_Q_scaled



Results Log

Name: StartX, Value: 30, Default?:
Name: EndX, Value: 55, Default?:
Name: Normalise, Value: 0, Defau
Name: Exclude, Value: , Default?:

Fit status: success
Stopped after 2 iterations
Fit successful, Duration 0.00 seconds

Algorithms

Execute Scale

Algorithms

- ▶ Arithmetic
- ▶ CorrectionFunctions
- ▶ Crystal
- ▶ DataHandling
- ▶ Diagnostics
- ▶ Diffraction
- ▶ Events
- ▶ Examples
- ▶ ILL
- ▶ Inelastic
- ▶ MDAlgorithms
- ▶ Muon

Workspaces

Load Delete Group

Filter Workspaces

Workspaces

- ▶ BaTiO3
- ▶ BaTiO3 Q
- ▶ BaTiO3 Q scaled


Script Interpreter

Fit Function (Chi-sq = 1.2008, failed)

Fit Display Setup

Property	Value
▼ Functions	
Type	CompositeFunction
NumDeriv	☐ False
▼ f0.FlatBackground	
Type	FlatBackground
A0	0.683681
▼ Settings	
Workspace	BaTiO3 Q scaled
Workspace Index	0
StartX	30.000000
EndX	55.000000
Output	BaTiO3 Q scaled
Minimizer	Levenberg-Marquardt
Ignore invalid ...	☐ False
Cost function	Least squares
Max Iterations	500
Peak Radius	0
Plot Difference	☒ True
Plot Composite ...	☐ False
Convolve Comp...	☐ False
Show Paramete...	☐ False
Evaluate Functi...	CentrePoint

BaTiO3_Q_scaled



Status: cannot reach the specified tolerance in F

Script Interpreter

%quickref -> Quick reference.

Results Log

Name: StartX, Value: 30, Default?: No, Direction: Input
Name: EndX, Value: 55, Default?: No, Direction: Input
Name: Normalise, Value: 0, Default?: Yes, Direction: Input
Name: Exclude, Value: , Default?: Yes, Direction: Input

Fit status: cannot reach the specified tolerance in F

Stopped after 6 iterations

Fit successful, Duration 0.00 seconds

Algorithms

Execute Scale

Algorithms

- ▶ Arithmetic
- ▶ CorrectionFunctions
- ▶ Crystal
- ▶ DataHandling
- ▶ Diagnostics
- ▶ Diffraction
- ▶ Events
- ▶ Examples
- ▶ ILL
- ▶ Inelastic
- ▶ MDAlgorithms
- ▶ Muon

Workspaces

- Load Delete Group Sort Save
- Filter Workspaces
- Workspaces
- ▶ BaTiO3
 - ▶ BaTiO3 Q
 - ▶ BaTiO3 Q scaled
 - ▶ BaTiO3 Q scaled NormalisedCovarianceMatrix
 - ▶ BaTiO3 Q scaled Parameters
 - ▶ BaTiO3 Q scaled Workspace

play coordinates!

DejaVu Sans


10 B It U x² x₂ αβ Γ f

Fit Function (Chi-sq = 1.2008, failed)

Fit Display Setup

Property	Value
▼ Functions	
Type	CompositeFunction
NumDeriv	☐ False
▼ f0-FlatBackground	
Type	FlatBackground
A0	0.683681
▼ Settings	
Workspace	BaTiO3 Q scaled
Workspace Index	0
StartX	30.000000
EndX	55.000000
Output	BaTiO3 Q scaled
Minimizer	Levenberg-Marquardt
Ignore invalid ...	☐ False
Cost function	Least squares
Max Iterations	500
Peak Radius	0
Plot Difference	☒ True
Plot Composite ...	☐ False
Convolve Comp...	☐ False
Show Paramete...	☐ False
Evaluate Functi...	CentrePoint

BaTiO3-Q-scaled-1



Scale input dialog

Scales an input workspace by the given factor, which can be either multiplicative or additive.

InputWorkspace: BaTiO3_Q_scaled_to_1


OutputWorkspace: BaTiO3_Q_scaled_to_1

Factor: 1.4627

Operation: Multiply

Keep Open ☐ Run Close

BaTiO3-Q-scaled-to-1-1



Results Log

Name: InputWorkspace, Value: BaTiO3_Q_scaled, De
Input
Name: OutputWorkspace, Value: BaTiO3_Q_scaled_to
Direction: Output
Name: Factor, Value: 1.4626999999999999, Default?:
Name: Operation, Value: Multiply, Default?: Yes, Dire
Scale successful, Duration 0.00 seconds

Algorithms

Execute Scale

Algorithms

- ▶ Arithmetic
- ▶ CorrectionFunctions
- ▶ Crystal
- ▶ DataHandling
- ▶ Diagnostics
- ▶ Diffraction
- ▶ Events
- ▶ Examples
- ▶ ILL
- ▶ Inelastic
- ▶ MDAlgorithms
- ▶ Muon

Workspaces

Load Delete Group Sort Save

Filter Workspaces

Workspaces

- ▶ BaTiO3
- ▶ BaTiO3 Q
- ▶ BaTiO3 Q scaled
- ▶ BaTiO3 Q scaled NormalisedCovarianceMatrix
- ▶ BaTiO3 Q scaled Parameters
- ▶ BaTiO3 Q scaled to 1
- ▶ BaTiO3 Q scaled Workspace

Status: cannot reach the specified tolerance in F

Script Interpreter

%quickref -> Quick reference.
help -> Python's own help system.
object? -> Details about 'object', use 'object??' for extra details.
%uiref -> A brief reference about the graphical user interface

play coordinates!

DejaVu Sans

10

B

$\int$

U

x^2

x_2

$\alpha\beta$

Γ

f

Fit Function (Chi-sq = 1.2008, failed)

Fit

Display

Setup

Property	Value
▼ Functions	
Type	CompositeFunction
NumDeriv	☐ False
▼ f0.FlatBackground	
Type	FlatBackground
A0	0.683681
▼ Settings	
Workspace	BaTiO3_Q_scaled
Workspace Index	0
StartX	30.000000
EndX	55.000000
Output	BaTiO3_Q_scaled
Minimizer	Levenberg-Marquardt
Ignore invalid ...	☐ False
Cost function	Least squares
Max iterations	500
Peak Radius	0
Plot Difference	☒ True
Plot Composite ...	☐ False
Convolve Comp...	☐ False
Show Paramete...	☐ False
Evaluate Functi...	CentrePoint

Status: cannot reach the specified tolerance in F

Script Interpreter

%quickref -> Quick reference.

help -> Python's own help system.

BaTiO3-Q-scaled-1

PDFFourierTransform input dialog

Fourier transform from S(Q) to G(r), which is paired distribution function (PDF). G(r) will be stored in another named workspace.

InputWorkspace

BaTiO3_Q_scaled_to_1

OutputWorkspace

BaTiO3_Gr

Reciprocal Space

InputSofQType

S(Q)

Qmin

0.9

Qmax

35

☐ Filter

Real Space

PDFType

G(r)

DeltaR

0.02

Rmax

50

rho0

0.002

Keep Open

☒

Run

Close

1

Results Log

Name: InputWorkspace, Value: BaTiO3_Q_scaled, Default?: No, Di

Input

Name: OutputWorkspace, Value: BaTiO3_Q_scaled_to_1, Default?:

Direction: Output

Name: Factor, Value: 1.4626999999999999, Default?: No, Direction

Name: Operation, Value: Multiply, Default?: Yes, Direction: Input

Scale successful, Duration 0.00 seconds

Algorithms

Execute

PDFFourierTransform

Algorithms

▶ Arithmetic

▶ CorrectionFunctions

▶ Crystal

▶ DataHandling

▶ Diagnostics

▶ Diffraction

▶ Events

▶ Examples

▶ ILL

▶ Inelastic

▶ MDAlgorithms

▶ Muon

Workspaces

Load

▼

Delete

Group

Sort

▼

Save

▼

Filter Workspaces

Workspaces

▶ BaTiO3

▶ BaTiO3 Q

▶ BaTiO3 Q scaled

▶ BaTiO3 Q scaled NormalisedCovarianceMatrix

▶ BaTiO3 Q scaled Parameters

▶ BaTiO3 Q scaled to 1

▶ BaTiO3 Q scaled Workspace

play coordinates!

DejaVu Sans

10

B

It

U

x²

x₂

αβ

Γ

f

Fit Function (Chi-sq = 1.2008, failed)

Ba

BaTiO3-Gr-1

PDFFourierTransform input dialog

Fourier transform from S(Q) to G(r), which is paired distribution function (PDF). G(r) will be stored in another named workspace.

InputWorkspace BaTiO3_Q_scaled_to_1

OutputWorkspace BaTiO3_Gr

Reciprocal Space

InputSofQType S(Q)

Qmin 0.9

Qmax 35

☐ Filter

Real Space

PDFType G(r)

DeltaR 0.02

Rmax 50



rho0 0.002

Keep Open ☒

Run

Close

BaTiO3_Gr



Status: cannot reach the specified tolerance in F

Script Interpreter

%quickref -> Quick reference.

help -> Python's own help system.

object? -> Details about 'object', use 'object??' for extra details.

Results Log

Name: Rmax, Value: 50, Default?: No, Direction: Input
Name: rho0, Value: 0.002, Default?: No, Direction: Input
Subtracting one from all values
Multiplying all values by Q
Adjusting to data: Qmin = 0.900289 Qmax = 35.0192
Using rmin = 0.02Angstroms and rmax = 50.02Angstroms
PDFFourierTransform successful, Duration 0.88 seconds

Algorithms

Execute PDFFourierTransform

Algorithms

- Arithmetic
- CorrectionFunctions
- Crystal
- DataHandling
- Diagnostics
- Diffraction
- Events
- Examples
- ILL
- Inelastic
- MDAlgorithms
- Muon

Workspaces

- Load Delete Group Sort Save
- Filter Workspaces
- Workspaces
- BaTiO3
 - BaTiO3_Gr
 - BaTiO3_Q
 - BaTiO3_Q_scaled
 - BaTiO3_Q_scaled NormalisedCovarianceMatrix
 - BaTiO3_Q_scaled Parameters
 - BaTiO3_Q_scaled to 1
 - BaTiO3_Q_scaled Workspace

MantidPlot - /SNS/groups/pg3-cycles/2017-A-B-data/PDF/Ni.mantid

File Edit View Graph Data Format Windows Catalog Interfaces Help

play coordinates! DeJaVu Sans 10 B It U x² x₂ αβ Γ f

Fit Function (Chi-sq = 1.2008, failed)

PDFFourierTransform input dialog

Fourier transform from S(Q) to G(r), which is paired distribution function (PDF). G(r) will be stored in another named workspace.

InputWorkspace: BaTiO3_Q_scaled_to_1

OutputWorkspace: BaTiO3_Gr

Reciprocal Space

InputSofQType: S(Q)

Qmin: 0.9

Qmax: 35

☐ Filter

Real Space

PDFType: G(r)

DeltaR: 0.02

Rmax: 50

rho0: 0.002

Keep Open ☒ Run Close

BaTiO3_Gr

SavePDFGui input dialog

Save files readable by PDFGui

InputWorkspace: BaTiO3_Gr

Filename: Browse

Keep Open ☐ Run Close

Save file

Name: BaTiO3.gr

Save in folder: SNS groups pg3-cycles 2017-A-B-data PDF Create Folder

Places: Search Recently Used 3ah Desktop File System Documents Downloads NOM PG3 CALIBRATION

Name: Ni Size: 10:34 Modified: 10:34

Cancel Save

Results Log

Name: Rmax, Value: 50, Default?: No, Direction: Input

Name: rho0, Value: 0.002, Default?: No, Direction: Input

Subtracting one from all values

Multiplying all values by Q

Adjusting to data: Qmin = 0.900289 Qmax = 35.0192

Using rmin = 0.02Angstroms and rmax = 50.02Angstroms

PDFFourierTransform successful, Duration 0.88 seconds

Algorithms

Execute SavePDFGui

Algorithms

- Arithmetic
- CorrectionFunctions
- Crystal
- DataHandling
- Diagnostics
- Diffraction
- Events
- Examples
- ILL
- Inelastic
- MDAlgorithms
- Muon

Workspaces

ad Delete Group Sort Save

Workspaces

- BaTiO3
- BaTiO3 Gr
- BaTiO3 Q
- BaTiO3 Q scaled
- BaTiO3 Q scaled NormalisedCovarianceMatrix
- BaTiO3 Q scaled Parameters
- BaTiO3 Q scaled to 1
- BaTiO3 Q scaled Workspace

Status: cannot reach the specified tolerance in F

Script Interpreter

- %quickref -> Quick reference.
- help -> Python's own help system.
- object? -> Details about 'object', use 'object??' for extra details.
- %gui ref -> A brief reference about the graphical user interface.