

## QSG – HYSPEC Crystal Alignment Guide

### Cheat sheet for aligning a crystal using a monochromatic neutron beam

1. If you're starting from scratch at HYSPEC, you're already too late!
  - a. CG1 and x-ray Laue alignment stations available.
    - i. X-ray Laue limitation: surface sensitive (vs. neutrons bulk sensitive)
    - ii. CG1 at HFIR is only available when reactor is on
  - b. Mounting sample
    - i. Wrap in Al foil before applying any wax. Wax is hard to remove and scatters many neutrons.
    - ii. We have a set of alignment jigs and tools that can be helpful aligning a single crystal.
    - iii. Know dimensions and mounting interfaces for your sample environment
    - iv. Mark faces with shiny sharpie pens
2. Anticipating Bragg reflections
  - a. Know basics of crystal before you arrive
  - b. Lower incident energy  $E_i$  restricts Bragg peaks to low HKL indices, reducing confusion
  - c. Location in scattering angle
    - i. ICSD website
    - ii. HYSPEC has some planning tools (Crystal class and WhereIsBragg)
    - iii. May need to move detector vessel to see desired scatter angle
  - d. Expected intensity
    - i. ICSD
    - ii. Others are available.
3. Finding the Bragg peak
  - a. Slow down sample rotation and start a run
  - b. Hopefully you can eyeball the right angle to within a few degrees to start
    - i. If needed, rotate  $\sim 180$  degrees and stop when you observe desired peak
    - ii. If not visible, tilt sample by  $\sim 6$  deg on one tilt axis and repeat, because the Bragg peak may be way out of plane
    - iii. Still not visible? Tilt even more or with other tilt axis
  - c. Once you are close, set a Region of Interest, and scan rotation angle  $\theta$ . Plot is generated with a  $\sim$ Gaussian peak. Center of peak is where desired Bragg peak is.
    - i. If the peak count rate is over  $\sim 300$  c/s, you run the risk of saturating a detector tube. You may need an attenuator.

- d. Confirm another peak about 90 degrees away.
- 4. Orienting a crystal plane in the horizontal plane
  - a. Set 'mirror' plane on tilt axis to avoid crosstalk
  - b. Cubic crystals: easy to find 90 deg separation between Bragg peaks (e.g. (100) and (010)); good for orthogonal tilt axes
  - c. Other crystals: look for roughly 90 deg Bragg peaks
- 5. Horizontal translation
  - a. Not needed if sample environment well mounted and sample well mounted
  - b. Increases flux.
  - c. Like tilt, translation x and y are better done on orthogonal reflections
- 6. Reducing background
  - a. Inelastic scattering:  $\sim 10^{-4}$  of elastic scattering, so background can be very bad
  - b. ToF background: elastic scatter away from the sample can look like energy gain or energy loss
  - c. Tighten up the upstream apertures using scans!