


TAX

Triple-Axis Spectrometer

HB-3 is a high-flux thermal neutron three-axis spectrometer designed for inelastic measurements on single crystals over a wide range of energy and momentum transfers. The energy and momentum range for measurements is quite large at HB-3. Because of its location directly at the end of the beam tube and the availability of a beryllium monochromator, the instrument is the ideal location for performing experiments at high-energy transfers (up to about 150 meV). The HB-3 monochromator provides three crystal choices (PG 002, Be 002, and Si 111) with variable vertical focus. This focus is calibrated to maintain a 2.54 cm beam height at the sample position, thus optimizing incident neutron flux as the incident energy varies. Of the three monochromators, pyrolytic graphite provides the highest neutron intensity as a result of its very high neutron reflectivity. The high-quality beryllium monochromator allows measurements with good energy resolution at higher energy transfers, whereas the silicon monochromator has the advantage of an absent second-order reflection, providing a higher order contamination-free beam.



SPECIFICATIONS

Beam spectrum	Thermal
Monochromators	Variable vertical focusing PG (002), Be (002), Si (111)
Analyzer	Fixed vertically focused PG (002)
Monochromator takeoff angle	$2\Theta_M = 12$ to 67°
Sample angle	$\pm 180^\circ$
Scattering angle	-90 to 128°
Analyzer angle	-40 to 90°
Detector	Single ^3He gas counter
Resolution (elastic)	5-10% E_i (adjustable with collimators)
Collimations (FWHM)	Premonochromator: 30', 48' Monochromator-sample: 20', 40', 60', 80' Sample-analyzer: 20', 40', 60', 80' Analyzer-detector: 30', 70', 90', 120', 210', 240'

26-G023702/cpp April 2026

APPLICATIONS

- Spin and lattice dynamics in high-temperature superconductors
- Magnetic excitations in low-dimensional magnetic model systems
- Spin waves in magnetically ordered systems
- Phonon anharmonicities and linewidths in thermoelectric materials and alloys
- Magnetic critical scattering in phase transitions
- Crystalline electric fields

For more information, contact

Songxue Chi, chis@ornl.gov, 865.603.3241Eleanor Clements, clementsem@ornl.gov, 865.863.5536neutrons.ornl.gov/taxManaged by UT-Battelle LLC
for the US Department of Energy