

Single Crystal Diffraction: The Definitive Structural Technique

Christine M. Beavers


18th National School on Neutron & X-ray Scattering

August 2nd, 2016

ALS

crys·tal·log·ra·phy

the branch of science dealing with the formation and properties of crystals




What is a Crystal?

- A crystal is a periodic arrangement of matter




What is a Crystal?

- A Crystal is a three-dimensional repeating array of atoms or molecules.
- In this example, our molecule is going to be in a shoebox, for simplicity.





Dimensionality





A one-dimensional array of shoeboxes

Dimensionality



A two-dimensional array of shoeboxes

Dimensionality





A three-dimensional array of shoeboxes

From Shoeboxes to Unit Cells



Unit Cell Definition:

Smallest volume unit of highest symmetry which, when translated in 3D will generate the crystal




The dimensions of the Unit Cell are an identifying feature for a specific crystal!

Crystal Selection #LifeGoals




Nice crystals are more likely to have nice diffraction


Crystal to Structure




Diffractometer Schematic




SXD Diffractometer





Data Collection Schematic



Crystal to Structure




Indexing





a=12.31Å, $\alpha=90.00^\circ$, V=1863Å³
b=12.31Å, $\beta=90.00^\circ$, Cubic I
c=12.31Å, $\gamma=90.00^\circ$

<input checked="" type="checkbox"/> Unit cell	
a [Å]	12.3052 ± 0.0008
b [Å]	12.3052
c [Å]	12.3052
α [°]	90.00
β [°]	90.00
γ [°]	90.00
V [Å ³]	1863.2 ± 0.4
<input type="checkbox"/> Domain translation	
x [mm]	0.04
y [mm]	-0.02


Crystal to Structure



Data Integration





Crystal to Structure




Absorption Correction



Overall scale and R(int) variations for Test



Crystal to Structure



Space Group Determination



Space Group Determination


3.1. SPACE-GROUP DETERMINATION AND DIFFRACTION SYMBOLS

Table 3.1.4.1. Reflection conditions, diffraction symbols and possible space groups (cont.)

ORTHORHOMBIC, Laue class mmm ($2/m\ 2/m\ 2/m$) (cont.)


Reflection conditions								Laue class mmm ($2/m\ 2/m\ 2/m$)			
hkl	$Ok\bar{l}$	$h0l$	$hk0$	$h00$	$Ok0$	$00l$	Extinction symbol	Point group	$mm2$ $m2m$ $2mm$	mmm	
	$k+l$	l	k		k	l	$Pncb$	222	$Pnn2$ (34)	$Pncb$ (50)	
	$k+l$	l	$h+k$	h	k	l	$Pncn$			$Pncn$ (52)	
	$k+l$	$h+l$		h	k	l	$Pnn-$			$Pnnm$ (58)	
	$k+l$	$h+l$	h	h	k	l	$Pnna$			$Pnna$ (52)	
	$k+l$	$h+l$	k	h	k	l	$Pnnb$			$Pnnb$ (52)	
	$k+l$	$h+l$	$h+k$	h	k	l	$Pnnn$			$Pnnn$ (48)	
$h+k$	k	h	$h+k$	h	k		$C---$			$C222$ (21)	$Cmm2$ (35) $Cm2m$ (38) $C2mm$ (38)
$h+k$	k	h	$h+k$	h	k	l	$C-2_1$	$C222_1$ (20)	$Cm2e$ (39)	$Cmme$ (67)	
$h+k$	k	h	h, k	h	k		$C-(ab)$				$C2me$ (39)
$h+k$	k	h, l	$h+k$	h	k	l	$C-c-$	$C222_1$ (20)	$Cmc2_1$ (36)	$Cmcm$ (63)	
$h+k$	k	h, l	h, k	h	k	l	$C-c(ab)$				$C2cm$ (40)
$h+k$	k	h, l	h, k	h	k	l					$C2ce$ (41)

Crystal to Structure



Electron Density from Diffraction

$$\rho_{xyz} = \frac{1}{V} \sum_{hkl} \bar{F}_{hkl} e^{-i2\pi(hx+ky+lz)}$$



Structure Solution


$$\rho_{xyz} = \frac{1}{V} \sum_{hkl} \bar{F}_{hkl} e^{-i2\pi(hx+ky+lz)}$$

-5	10	-3	0.01	0.45
-6	10	-3	-0.28	0.49
-7	10	-3	-0.28	0.52
-8	10	-3	1.63	0.56
12	9	-3	1.15	0.68
-12	-9	3	0.64	0.66
11	9	-3	5.65	0.82
-11	-9	3	6.17	0.77
10	9	-3	-0.14	0.65
-10	-9	3	-0.14	0.65
9	9	-3	2.01	0.72
-9	-9	3	2.11	0.51
8	9	-3	4.10	0.72
8	9	-3	6.25	0.82
-8	-9	3	5.89	0.63
7	9	-3	9.05	0.95
-7	-9	3	9.79	0.79
6	9	-3	4.40	0.72
-6	-9	3	7.25	0.65
5	9	-3	8.03	0.84
-5	-9	3	7.17	0.73
4	9	-3	5.66	0.71
-4	-9	3	4.98	0.65
3	9	-3	1.28	0.51
-3	-9	3	1.40	0.46
2	9	-3	18.66	1.23
-2	-9	3	16.45	1.18
1	9	-3	8.06	0.84
-1	-9	3	7.75	0.75
0	0	-3	17.88	1.22


$$F_{hkl} = |F_{hkl}| e^{i\phi_{hkl}}$$

????


SHELXT
SHELXS
SIR2011+
SUPERFLIP





What is a Structure?



Crystal to Structure



Refinement & Validation





$$R_1 = \frac{\sum_{hkl} ||F_{obs}| - |F_{calc}||}{\sum_{hkl} ||F_{obs}||}$$

$$wR_2 = \sqrt{\frac{\sum_{hkl} w(F_{obs}^2 - F_{calc}^2)^2}{\sum_{hkl} w(F_{obs}^2)^2}}$$

More Resources!!!

Books

- Werner Massa
 - ISBN-13: 978-3-540206446
- Stout & Jensen
 - ISBN-13: 978-0471607113



Internet

- X-ray Forum
 - www.xrayforum.co.uk/
- IUCr Forum
 - forums.iucr.org
- CCP4
 - <http://www.ccp4.ac.uk>


Small Molecule Crystallography at a Synchrotron *or* What can you do with more flux?

Contents

- Why do crystals diffract poorly?
- What can we do to them to make them diffract poorly?
- What can we learn from poorly diffracting crystals?
- What do synchrotrons have to do with all this?


The Spectrum of Crystallinity

- Perfect crystals
 - Diffract dynamically, Bragg



The Spectrum of Crystallinity


- Good Crystals
 - Diffract kinematically(Bragg), due to mosaicity, but still have good long range order




The Spectrum of Crystallinity

- Poor crystals

- Diffract kinematically(Bragg), but diffraction limited due to poor long range order.
- Can show powder Laue rings/spot smearing due to mosaicity becoming microcrystallinity
- Can also display non-Bragg scatter due to TDS



Positional Disorder



Poor Long Range Order


Scattering Efficiency

$$\text{Intensity of Diffraction} \approx \lambda^3 = \frac{LI_{\text{incident}} \langle |F_{hkl}^2| \rangle V_{\text{crystal}}}{V_{\text{cell}}^2}$$


- where:-
- F = number of electrons per atom
- V_{crystal} = volume of the crystal
- V_{cell} = volume of the unit cell

M.M.Harding *J. Synchrotron Radiation*, 250-259 1996


Effect of disorder



Intensity vs. Displacement



Higher Angle Reflections Affected by Larger ADPs



Wavelength

- The material and the wavelength need to be compatible
 - Short wavelengths better for heavy absorbers
 - Long wavelengths better for light atoms (weakly diffracting elements)
 - Be aware of absorption edges and potential fluorescence from sample

Bigger isn't always better

- Large crystals aren't guaranteed to diffract better
- Crystal should match beam size
 - But if there is a choice, smaller than the beam is usually better
- Rocking width can be worse with large crystals due to poor mosaicity


Structures from change:

IN SITU EXPERIMENTS

In-situ Crystallography

- The application of a stimuli to produce structural change
 - Temperature
 - Pressure
 - Gas or Vacuum
 - Light
 - Electric or Magnetic Fields

Desolvation




Three-Way Crystal-to-Crystal Reversible Transformation and Controlled Spin Switching by a Nonporous Molecular Material

Sanchez Costa et al., *J. Am. Chem. Soc.*, **2014**, 136 (10), pp 3869–3874


DOI: 10.1021/ja411595y

ALS

Photocrystallography



Experimental Procedure



Procedure

- High quality ground state data collection
 - Irradiation (LEDs) - LED ring
- Metastable state data collection
 - Inspection of the density map
- Temperature variation experiments



N-bound

O-bound

nitro

nitrito

0 %


50 %

100 %

No irradiation



10 minutes irradiation

20 minutes irradiation




J. Appl. Cryst. 2010, **43**, 337-340


Gas Cell



Hydrated MOF




Dehydrated




ALS


NO absorbed



SO₂ absorbed



High Pressure with Diamond Anvil Cells



Diamond Anvil Cell (DAC)





Figure 2
Exploded view of a Merrill-Bassett DAC.

Moggach, S. A. et al. *J. Appl. Crystallogr.* **2008**, *41*, 249-251.

ALS

A High Pressure Sample



The Zinc-Alkyl Gate (ZAG) Family





Figure 4. Comparison of ZAG-4 (left) and ZAG-6 (right) as viewed down their c -axes.

ZAG-4, under pressure



Why Synchrotrons?

- In situ experiments usually produce the degradation of a crystal, and most are more successful with small crystals.
- Poorly diffracting crystals need as much intensity as they can take.
- In both cases, a synchrotron offers orders of magnitude more flux, which means a better chance of success.

ALS

The Advanced Light Source is supported by the Director, Office of Science, Office of Basic Energy Sciences, of the U.S. Department of Energy under Contract No. DE-AC02-05CH11231

COMPRES, the Consortium for Material Properties Research in Earth Sciences, is supported under NSF Cooperative Agreement EAR 11-57758.

Advanced Light Source
An Office of Science User Facility

Thank You!!

